



## *Washing Buffer based on TRIS without detergents* **For use in immunohistochemistry**

### Product description

Composition:	TRIS-based buffer with 350 mM NaCl 10x stock solution contains neither Tween nor other detergents
pH-value:	pH 7.2 ± 0.2
Preservative:	contains 0.1% ProClin® 300
Storage:	10x stock solution: shelf life at –20 °C: 1 year - tolerates repeated freezing and thawing cycles shelf life at 2 – 8 °C: 6 months

Use working solution immediately!

**For research use only, not for diagnostic use**

### Available

Package sizes: 500 ml order number 146 500

### Instructions for use

This *Washing Buffer* is designed especially for immunohistochemistry. It should be used, whenever detergents can make problems in detection.

Crystals of salt can precipitate after storage at 2 – 8 °C or after freezing. Therefore *Washing Buffer* must be warmed up to room temperature and should be mixed thoroughly before preparing the working solution. This leads to dissolving of salt after shaking. Stock solution is diluted 1:10 with salt free water to get the working solution.

Use working dilution immediately.

For further information please visit [www.candor-bioscience.com](http://www.candor-bioscience.com).

ProClin is a registered trade mark of Rohm and Haas Company.