



## *Antibody Stabilizer*

### **Stabilizer for long-term storage of proteins or antibodies at 2-8°C**

Available products: *Antibody Stabilizer TRIS* (article no. 130)  
*Antibody Stabilizer PBS* (article no. 131)

Storage: 2-8°C (Does not tolerate freezing!)

pH-value: 7.2 ± 0.2 (article no. 130)  
7.4 ± 0.2 (article no. 131)

Preservative: contains < 0.0014 % [w/w] reaction mass of CMIT/MIT (3:1)

Expiry date  
when stored unopened: please refer to the label on the bottle

**For research use only, not for diagnostic use**

### **Instructions for use**

*Antibody Stabilizer* is ready-to-use. Please shake the buffer thoroughly before use.

Just dissolve proteins or antibodies with *Antibody Stabilizer* for storage in the refrigerator.

Antibodies/proteins can be diluted 1:20 or in lower concentrations for storage.  
Typical concentrations for stored antibodies are between 400 and 1000 ng/mL.  
A lot of antibodies can also be stored for long-term in *Antibody Stabilizer* without loss of binding activity in very low concentrations (e.g. 80 ng/mL)  
Low concentrations during storage save pre-dilution steps before use of the antibodies.

Stability data of one antibody or protein cannot be transferred to other antibodies/proteins. Therefore every antibody/protein has to be tested for its shelf life in *Antibody Stabilizer*. If *Antibody Stabilizer* is used for immunodiagnostic kits, the shelf life has to be tested according to the regulatory requirements for diagnostics, which apply.

Please note that high protein concentration and/or microbial contamination may lead to reduced effectivity of the preservative. If you add your conjugate/antibody for storage and you are not sure about the potential microbial contamination, it may be beneficial to add additional preservatives or antibiotics.

For further information please visit [www.candor-bioscience.com](http://www.candor-bioscience.com).