

- **Antibody Stabilizer**
- Liquid Plate Sealer®
- HRP-Protector™

Antibody Stabilizer

Stabilizer for long-term storage of proteins and antibodies in solution. Trouble-free long-term storage at 4°C!

Antibody Stabilizer contains additives for preserving the conformation of proteins and antibodies in solution.

Antibody Stabilizer

- stabilises by conserving protein structure
- optionally based on TRIS or PBS



Antibody Stabilizer conserves the structure of proteins and antibodies preventing them from losing functionality during storage. Just dissolve proteins or antibodies in Antibody Stabilizer and store in the refrigerator.

Ready-to-use


Available package sizes:

Antibody Stabilizer based on TRIS

50 ml order number 130 050
 125 ml order number 130 125
 500 ml order number 130 500

Antibody Stabilizer based on PBS

50 ml order number 131 050
 125 ml order number 131 125
 500 ml order number 131 500



Shelf life of an antibody is measured as binding in an ELISA after storage under temperature stress. The tests were performed with an antibody-concentration of only 400 ng/ml.