

- The Blocking Solution
- SmartBlock™
- **BSA-Block**



BSA-Block

Blocker for ELISA, EIA, Western blotting, protein arrays and immuno-PCR.

BSA-Block

- saturates free binding sites
- well established
- competitively priced



BSA-Block saturates free binding capacities on plastic consumables and other surfaces like ELISA plates and blotting membranes.

BSA-Block is a well established, versatile blocker for many applications.

The high purity and control of the raw materials enable good lot-to-lot consistency for this blocker based on BSA. Therefore a high reliability of the assay over years is given even when different lots are used.

If you experience background due to non-specific binding even when using BSA-Block, we recommend our much more effective product The Blocking Solution.

Ready-to-use

Available package sizes:

125 ml order number 115 125

500 ml order number 115 500