

- Antibody Stabilizer
- Liquid Plate Sealer®
- **HRP-Protector™**

HRP-Protector™

HRP-Protector™ preserves enzymatic activity and conformation of your HRP (horse radish peroxidase) conjugates.

- for long-term stability
- ready-to-use



HRP-Protector™ gives peroxidase conjugates the crucial long-term stability over several years even when stored in low concentrated ready-to-use dilutions. It maintains the activity of conjugates of antibodies and of Neutravidin/Streptavidin with HRP.

HRP-Protector™ can be used directly as an assay buffer for incubating the HRP-conjugates. This enables assays without pre-dilution steps for the conjugates after long-term storage. No need to freeze and thaw, no need for pre-dilution.

Conjugates are simply stored ready-to-use at 4°C for easy handling.

Shelf life of HRP-antibody conjugate is measured as peroxidase signal in an ELISA after binding the analyte. Values are shown as percentage of freshly diluted HRP-antibody conjugate.

Ready-to-use

Available package sizes:

- 50 ml order number 220 050
- 125 ml order number 220 125
- 500 ml order number 220 500

Heat stress test at 45°C

