

- LowCross-Buffer®
- LowCross-HRP®
- The Blocking Solution

LowCross-HRP®

Gives you two excellent products in one great buffer solution.

- HRP enzyme stabilizer
- HAMA-blocker substitute
- Interference reducer



LowCross stands for superior reliability by reducing assay interferences including HAMA-derived effects.


LowCross-HRP® gives peroxidase-antibody conjugates the crucial long-term stability over several years even in a ready-to-use dilution.

Ready-to-use

Available package sizes:

- 50 ml order number 200 050
- 125 ml order number 200 125
- 500 ml order number 200 500

Heat stress test at 45°C



Shelf life of HRP-antibody conjugate is measured as peroxidase signal in an ELISA after binding the analyte. Values are shown as percentage of freshly diluted HRP-antibody conjugate.