

- **Sample Buffer**
- Washing Buffer
- Coating Buffer
- Stripping Buffer



# Sample Buffer

Sample and antibody dilution buffer for ELISA, EIA, RIA, protein arrays, Western blotting, immuno-PCR and immunohistochemistry.

- for diluting the analyte or specimen
- for diluting the detection antibody
- free of phosphate



Sample Buffer is the standard specimen and antibody diluent for many assays.

Sample Buffer is designed for trouble free assays without interference. If you have background or false-positives – for example derived from cross-reactivities or matrix effects - we recommend to use LowCross-Buffer® rather than Sample Buffer.

Specimen with analyte – as well as the detection antibody – is diluted in Sample Buffer and then used in the assay.

Sample Buffer is also available without Tween or other detergents for use in immunohistochemistry.

## Ready-to-use

Available package sizes:

### Sample Buffer

- 50 ml order number 105 050
- 125 ml order number 105 125
- 500 ml order number 105 500

### Sample Buffer

#### without Tween for IHC

- 50 ml order number 106 050
- 125 ml order number 106 125
- 500 ml order number 106 500