

- Sample Buffer
- Washing Buffer
- Coating Buffer
- **Stripping Buffer**

Stripping Buffer

For reprobing on Western blotting membranes.
Easy and safe removal of antibodies!

- ✓ convenient to use
- ✓ non-malodorant and non-hazardous,
because free from DTT and mercaptoethanol
- ✓ no need for a fume hood



Stripping Buffer removes reaction solution and primary and secondary antibodies from Western blotting membranes. After stripping the membrane can be used for repeated detection (reprobing) with different antibodies.

Stripping Buffer does not contain mercaptoethanol nor DTT.

Stripping Buffer is odorless and working in the lab becomes much more pleasant than with mercapto-ethanol based formulations!

Ready-to-use

Available package sizes:

50 ml order number 150 050
125 ml order number 150 125
500 ml order number 150 500