

- Sample Buffer
- **Washing Buffer**
- Coating Buffer
- Stripping Buffer



Washing Buffer

Wash buffer for ELISA, EIA, RIA, protein arrays, Western blotting, immuno-PCR and immunohistochemistry.

Washing Buffer

- ✓ time saving
- ✓ optionally based on TRIS or PBS
- ✓ suitable for ELISA washer

Be more dependable



Washing steps are needed to remove unbound and excessive components which are able to interfere with the assay. Washing Buffer is suitable for automatic washers. (Please note the specifications of your device and follow the instruction manual of your washer!)

Buffer is optionally based on TRIS or PBS.

For immunohistochemistry also available without Tween and detergents.

10 x stock solution
Available package sizes:

Washing Buffer PBS
500 ml order number 140 500

Washing Buffer PBS without Tween for IHC
500 ml order number 141 500

Washing Buffer TRIS
500 ml order number 145 500

Washing Buffer TRIS without Tween for IHC
500 ml order number 146 500